# Articles

# Oral swabs with a rapid molecular diagnostic test for pulmonary tuberculosis in adults and children: a systematic review

E Chandler Church, Karen R Steingart, Gerard A Cangelosi, Morten Ruhwald, Mikashmi Kohli, Adrienne E Shapiro

# **Summary**

**Background** Tuberculosis is a leading cause of infectious disease mortality worldwide, but diagnosis of pulmonary tuberculosis remains challenging. Oral swabs are a promising non-sputum alternative sample type for the diagnosis of pulmonary tuberculosis. We aimed to assess the diagnostic accuracy of oral swabs to detect pulmonary tuberculosis in adults and children and suggest research implications.

Methods In this systematic review, we searched published and preprint studies from Jan 1, 2000, to July 5, 2022, from eight databases (MEDLINE, Embase, Scopus, Science Citation Index, medRxiv, bioRxiv, Global Index Medicus, and Google Scholar). We included diagnostic accuracy studies including cross-sectional, cohort, and case-control studies in adults and children from which we could extract or derive sensitivity and specificity of oral swabs as a sample type for the diagnosis of pulmonary tuberculosis against a sputum microbiological (nucleic acid amplification test [NAAT] on sputum or culture) or composite reference standard.

Findings Of 550 reports identified by the search, we included 16 eligible reports (including 20 studies and 3083 participants) that reported diagnostic accuracy estimates on oral swabs for pulmonary tuberculosis. Sensitivity on oral swabs ranged from 36% (95% CI 26–48) to 91% (80–98) in adults and 5% (1–14) to 42% (23–63) in children. Across all studies, specificity ranged from 66% (95% CI 52–78) to 100% (97–100), with most studies reporting specificity of more than 90%. Meta-analysis was not performed because of sampling and testing heterogeneity.

Interpretation Sensitivity varies in both adults and children when diverse methods are used. Variability in sampling location, swab type, and type of NAAT used in accuracy studies limits comparison. Although data are suggestive that high accuracy is achievable using oral swabs with molecular testing, more research is needed to define optimal methods for using oral swabs as a specimen for tuberculosis detection. The current data suggest that tongue swabs and swab types that collect increased biomass might have increased sensitivity. We would recommend that future studies use these established methods to continue to refine sample processing to maximise sensitivity.

Funding Bill and Melinda Gates foundation (INV-045721) and FIND (Netherlands Enterprise Agency on behalf of the Minister for Foreign Trade and Development Cooperation [NL-GRNT05] and KfW Development Bank, German Federal Ministry of Education and Research [KFW-TBBU01/02]).

Copyright © 2023 The Author(s). Published by Elsevier Ltd. This is an Open Access article under the CC BY 4.0 license.

# Introduction

Tuberculosis remains a leading cause of infectious disease mortality and is among the top 15 causes of death worldwide.<sup>1</sup> An estimated 7.5 million people were diagnosed with tuberculosis in 2022, an increase from 2021. There were an estimated 1.3 million deaths from tuberculosis in 2022, with around 167000 of those deaths occurring in people living with HIV.<sup>1</sup> In 2022, only 63% of people treated for tuberculosis had a microbiological test confirming the diagnosis.<sup>1</sup> WHO recommends using a molecular rapid diagnostic test as the initial test for diagnosing tuberculosis in adults and children with signs and symptoms of tuberculosis.<sup>2</sup>

Specimens for diagnosing pulmonary tuberculosis include sputum, gastric aspirates, nasopharyngeal aspirates, urine, and stool. For diagnosing pulmonary tuberculosis using sputum, Xpert MTB/RIF and XPERT MTB/RIF Ultra represent a high sensitivity molecular testing method for adults and children when the specimen can be obtained.<sup>3,4</sup> Sputum production can be limited in individuals with a weak cough, people living with HIV, and children younger than 5 years. These limitations often require the use of other respiratory specimens (eg, gastric aspirates) in young children or induced sputum in adults, which require medical equipment and specialised staff training.<sup>5</sup>

To improve access to tuberculosis testing, oral swabs have been proposed as a new clinical specimen for pulmonary tuberculosis diagnosis. Oral swab specimens are analysed by nucleic acid amplification tests (NAAT) that amplify *Mycobacterium tuberculosis* DNA collected by the swab. Multiple types of swabs and NAAT platforms are being studied, with PCR the most commonly used NAAT method.<sup>5</sup> To date, there are no





#### Lancet Glob Health 2024; 12: e45–54

See Comment page e4 HIV Vaccine Trials Network, Fred Hutchinson Cancer Center, Seattle, WA, USA (E C Church MD); Division of Allergy and Infectious Diseases, University of Washington, Seattle, WA, USA (E C Church, A F Shapiro MD): Honorary Research Fellow, Department of Clinical Sciences, Liverpool School of Tropical Medicine, Liverpool, UK (K R Steingart MD); Department of Environmental and Occupational Health Sciences. School of Public Health, University of Washington, Seattle, WA, USA (Prof G A Cangelosi PhD): FIND. Geneva, Switzerland (M Ruhwald MD, M Kohli PhD); Department of Global Health, University of Washington. Seattle, WA, USA (A E Shapiro)

Correspondence to: Dr E Chandler Church, HIV Vaccine Trials Network, Fred Hutchinson Cancer Center, Seattle 98109, WA, USA cchurch2@fredhutch.org

#### **Research in context**

#### Evidence before this study

Despite recent advances in tuberculosis diagnostics, most methods still require sputum for testing. However, many populations, including children and people living with HIV, have difficulty producing sputum samples. Obtaining sputum for diagnosis in people who cannot produce it themselves can be resource intensive and require trained staff, which can make obtaining a microbiological diagnosis much more challenging. We conducted a search on July 5, 2022, in the following databases: Medline, Embase, Scopus, Science citation index, medRxiv, bioRxiv, Global index medicus, and Google Scholar using search terms including: (("Mycobacterium tuberculosis" OR "tuberculosis" OR "TB" OR "pulmonary tuberculosis") AND ("oral swab\*" OR "tongue swab\*" OR "buccal swab\*")), for systematic reviews looking at oral swabs. The search was repeated on Aug 8, 2023. We did not find any previous systematic reviews looking at oral swabs for tuberculosis at the time of the original search, and no reviews focusing specifically on oral swab methods are currently in existence, although Savage and colleagues addressed oral swab methods in a review of non-sputum-based diagnostics.

## Added value of this study

To our knowledge, this is the first comprehensive review of oral swabs as a specimen type for the diagnosis of pulmonary tuberculosis. The methods of the included studies were highly variable, particularly regarding the area of the oral cavity swabbed, the type of swab used, the processing of the swab sample before analysis, and the type of molecular testing performed on the swab. Several studies directly compared one of these factors. This review brings together all the available evidence on oral swabs to help future studies use the techniques that yield the best results in terms of sensitivity and specificity.

#### Implications of all the available evidence

Our study has substantial implications for future research into oral swabs. By bringing together the available evidence, researchers can choose the techniques that will maximise the potential sensitivity and specificity of oral swabs within the study population. We would recommend that researchers focus their efforts on tongue swabs and swab types that collect increased biomass, as these methods appear to increase sensitivity. In addition, although in-house PCR testing was primarily used in the early stages of research on oral swabs, recent studies show that reasonable sensitivity can be achieved using Xpert MTB/Rif Ultra, which is widely used for sputum testing and available to most sites that would benefit from this type of sample.

oral swab-based recommendations for tuberculosis detection by WHO.

Diagnosing pulmonary tuberculosis using oral swabs as a clinical specimen comprises a multi-stage process of sampling, processing, and analysis. Differences in each component of the analytic pathway can contribute to overall variability in diagnostic accuracy of oral swabs (figure 1). There are several potential advantages to this clinical specimen over sputum-based sampling, including possible self-collection and use in populations that can have difficulty producing sputum.<sup>6</sup> Oral swab samples could also be combined with other sample types such as Determine TB LAM Ag (Abbott, Lake Forest, IL, USA) or gastric aspirate in populations with low bacillary burden.<sup>7</sup> One potential clinical pathway for the use of oral swabs is shown in the appendix (p 8).

See Online for appendix

We performed a systematic review of existing oral swab literature to determine the diagnostic accuracy of this emerging sample type, identify sources of variability in accuracy, and characterise priorities for further optimisation and research.

## Methods

### Search strategy and selection criteria

We searched eight databases (MEDLINE, Embase, Scopus, Science Citation Index, medRxiv, bioRxiv, Global Index Medicus, and Google Scholar) for published and preprint studies from Jan 1, 2000, to July 5, 2022, without language restriction (see appendix p 1 for the full search strategy). The review was performed on July 5, 2022. Search strategy included terms for tuberculosis, oral swabs, and brand names for various oral swabs. We reviewed reference lists of included articles and any relevant review articles identified. We also contacted researchers at FIND and other experts in the field of tuberculosis diagnostics (including GAC and MK) for information on unpublished studies. Ethical approval was not required for this review.

We included cross-sectional, cohort, and case-control studies evaluating the sensitivity and specificity of oral swab sampling and molecular testing against a microbiological or composite reference standard for the diagnosis of pulmonary tuberculosis. The microbiological reference standard included solid or liquid culture or a molecular WHO-recommended rapid diagnostic performed on sputum. The composite reference standard combined the microbiological reference standard with clinical criteria including symptoms and chest imaging. We considered both published articles and preprints from which we could extract or derive true positive, false positive, false negative, and true negative values. We only included studies where both the reference standard and test strategy were assessed in all participants. We excluded studies that evaluated saliva as a testing sample as this review is focused on swab-based samples, which are not used for saliva collection. Oral swabs are collected using a porous or scraping material that is rubbed over a specific location in the oral cavity, such as the tongue or buccal mucosa, so that friction lifts surface cellular



Figure 1: Potential sources of variability in sensitivity and specificity for oral swabs

LAMP=loop-mediated isothermal amplification. Mfr=manufacturer. NAAT=nucleic acid amplification test.

material to the swab as well as absorbing oral fluids, whereas saliva samples are oral fluids deposited directly into a collection vessel, without use of an intermediary device. Participants were adults (aged 15 years or older) and children (younger than 15 years) who were evaluated for presumptive pulmonary tuberculosis.<sup>89</sup>

## Data analysis

We used Covidence systematic review software to manage the selection of studies.<sup>10</sup> Two review authors independently scrutinised titles and abstracts identified from literature searching to identify potentially eligible studies. We retrieved the article of any citation, identified by any review author, for full-text review. Then, two review authors (ECC and AES) independently assessed articles for inclusion using the predefined selection criteria. We resolved any discrepancies by discussion or with a third review author. We followed review methods described in the Cochrane Handbook for Systematic Reviews for Diagnostic Test Accuracy.<sup>11</sup>

Three review authors (ECC, KRS, and AES) independently extracted data on key participant characteristics and information on the index text (location of oral swabbing, type of swab, and type of NAAT), reference standard, and diagnostic  $2\times 2$  table data (number of true positives, false positives, false negatives, and true negatives). We used the QUADAS-2 tool, tailored to this review, to assess the quality of the included studies (appendix p 4).<sup>12</sup> Three review authors (ECC, KRS, and AES) independently assessed all domains for risk of bias and the first three domains for concerns regarding applicability. Study authors were contacted for clarifications and missing data. Review authors who were also authors of included studies did not perform quality assessment or extract data from their own study or studies.

We determined sensitivity and specificity estimates and 95% CIs for individual studies, and generated forest plots using Review Manager 5 (Review Manager 2020). We visually inspected forest plots for heterogeneity.

# Role of the funding source

The funding source had no role in the review of the literature, data abstraction, or quality assessment.

# Results

We obtained 965 records and, after deleting 415 duplicates, we had 550 records to screen. We excluded 509 records



#### Figure 2: Study selection

mWRD=molecular WHO-recommended rapid diagnostic.

based on title and abstract review and identified 41 for full-text review. We were unable to obtain a full-text copy of one study and were unable to contact the author, so this study was excluded.<sup>13</sup> We excluded 24 reports due to ineligibility. We identified 16 reports that met our inclusion criteria (figure 2). Four reports reported more than one oral swab testing strategy study, hence we included 20 studies (3083 participants) in the review. Participant characteristics for each study are summarised in the table. When different methods were evaluated in the same study (ie, different swab types), we used the designation a, b; for example, Flores 2020a and Flores 2020b used different swab types and were considered as two distinct studies.

Type of NAAT	Xpert MTB/Rif Ultra	Xpert MTB/Rif Ultra	In-house PCR	In-house PCR (Roche 480 Lightcycler)	In-house PCR (Qiagen TAQ PCR Core Kit)	In-house PCR (IS 6110-targeted qPCR [Applied Biosystems StepOnePlus Real- Time PCR system])	Xpert MTB/RIF Ultra Ultra inues on next page)
Laboratory processing	Method a: 2 swabs processed together with Cepheid sample reagent, then prepared for Xpert testing by the same protocol as sputum; methodb: 1 swab bied, incubated, and mixed without use of sample reagent use of sample reagent	Thawed and 1.6 mL Xpert sample reagent added, vortexed, then incubated	Vortexed with 2 mm glass beads	OmniSwabs vortexed before freezing, DNA extracted with QlAmp DNA minikit, EasiCollect cut into pieces and placed in pieces and placed in followed by neutral solution to release DNA	Samples liquified and injected into SLIM platform	Not specified	Not specified (Table cont
Sample storage	2 swabs processed within 1 h, remaining swab frozen and stored at -80°C until testing	Frozen and stored at -80°C until testing	Stored at 4°C	OmniSwab frozen and stored at -80°C, EasiCollect at -20°C	Frozen and stored at -80°C until testing	Frozen and stored at -80°C until testing	1 swab analysed within 24 h, second stored at -80°C for up to 12 weeks
Sample transport	2 swabs placed together in 800 mL of Tris-EDTA buffer. If additional swab collected, stored in a separate container with the same media type	Sterile saline	3 mL Middlebrook 7H9 with OADC and Tween 80	OmniSwab in sterile lysis buffer, EasiCollect dry	Not specified	Sterile lysis buffer	Eppendorf and sterile lysis buffer
Number, location of swabs	2–3 tongue swabs	1-2 buccal or tongue swabs	1 tongue swab	1 buccal swab	1 gingival swab	2 buccal swabs	2 tongue swabs 24 h apart
Swab type	Copan FLOQswab	Puritan Purflock and Copan FloQswab	Copan 480 C swab	Whatman OmniSwab or Whatman EasiCollect FTA card	OmniGene Oral OMR- 110 kit swab	Whatman OmniSwab	Puritan Purflock Ultra swabs
Positive by oralswab and reference standard	42 (Andama 2022a); 27 (Andama 2022b)	20	11	S	65	13	66
Reference standard (primary)	Sputum Xpert MTB/ RIF Ultra	Sputum Xpert MTB/ RIF Ultra or culture	Clinical diagnosis	Sputum or gastric aspirate culture	Sputum culture	Sputum Xpert MTB/ RIF or culture	Sputum Xpert MTB/ RIF
Participants living with HIV (%)	58 (31.7%)	57 (19-6%)	6 (17·1%)	o	1 (0.4%)	54 (54%)	Not reported
Children aged <5 years	:	Not reported	35 (100%)	124 (43%)	:	:	:
Median age (years)	33 (IQR 26-43)	2:7 (IQR 1:2-6:1)	0.4 (IQR 0.4- 1.6)	Not reported	Mean 58·8 (SD 15·2)	38 (IQR 30-44)	Not reported
Total partici- pants	184 (Andama 2022a); 37 (Andama 2022b)	291	35	288	272	100	256
Population	Adults aged >18 years presenting for health care with at least 2 weeks of cough	Children aged \$15 years with cough and one other finding suggestive of pulmonary tuberculosis	Children aged <5 years clinically diagnosed with tuberculosis or lower respiratory tract infection	Children aged <15 years with at least one symptom of pulmonary tuberculosis	Adults aged >18 years with presumptive tuberculosis	Adults aged ≥13 years with presumptive tuberculosis	Screening study in prison inmates, no age restrictions noted
	Andama 2022, Uganda <sup>22</sup>	Cox 2022, South Africa®	Ealand 2021, South Africa <sup>22</sup>	Flores 2020, Peru <sup>27</sup>	Kang 2022, South Korea <sup>20</sup>	LaCourse 2022, Kenya²³	Lima 2020, Brazil <sup>14</sup>

Type of NAAT		In-house PCR (IS 6110-targeted aPCR [Applied Biosystems StepOnePlus Real- Time PCR system])	Xpert MTB/RIF Ultra	In-house PCR (Roche 480 Lightcycler)	In-house PCR (Roche 480 Lightcycler)	In-house PCR	In-house PCR (IS 6110-targeted aPCR [Applied Biosystems StepOnePlus Real- Time PCR system]) system])	
Laboratory processing		DNA extracted with QIAmp DNA minikit	Not specified	FTA cards cut into pieces and placed in basic solution, then placed in neutral solution. QIAmp DNA minikit used to extract DNA. OmniSwabs warmed and extracted using QIAmp DNA minikit	DNA extracted with PrimeXtract nucleic acid extraction and purification kit	DNA extracted using QIAmp DNA minikit	Manually processed and concentrated with ethanol (Table cont	
Sample storage		Frozen and stored at -80°C within 8 h of collection	Unknown	OmniSwabs vortexed and stored at -80°C. FTA cards frozen at -20°C	Vortexed then kept at room temp for up to 30 days or frozen and stored at -20°C	Frozen and stored at -80°C	Frozen and stored at -80°C until testing	
Sample transport		Sterile lysis buffer	OmniSwab in sterile lysis buffer, EasiCollect dry	OmniSwab ejected into media, EasiCollect stored dry	PrimeStore molecular transport media	Sterile lysis buffer	0.5 mL stenie Tris- EDTA buffer	
Number, location of swabs		1 tongue swab with Purflock Ultra, 2 tongue swabs with Whatman OmniSwab	1 buccal swab	3 buccal swabs	1 buccal swab	2 buccal swabs	2 tongue swabs	
Swab type		Whatman OmniSwab and Purtan Purflock Ultra	Whatman OmniSwab or Whatman EasiCollect FTA card	Whatman OmniSwab or Whatman EasiCollect FTA card	PrimeSwab flocked swab	GE Healthcare OmniSwab or Puritan Purflock Ultra swab	Copan FLOQswab	
Positive by oralswab and reference standard		74 (Luabeya 2020a);74 (Luabeya 2020b) 2020b)	15	38 (Mesman 2020a); 15 (Mesman 2020b)	29	17	42	
Reference standard (primary)		Sputum Xpert MTB/ RIF	Sputum culture	Sputum culture	Sputum Xpert MTB/ RIF or culture	Sputum Xpert MTB/ RIF or culture or clinical diagnosis	Sputum Xpert MTB/ RIF Ultra or culture	
Participants living with HIV (%)		86 (25.1%)	Not reported	4 (3.3%)	0	18 (11%)	120 (92%)	
Children aged <5 years		:	AA	e z	AA	113(68.5%)	¥ Z	
Median age (years)		35-2 (QR 25-6- 45-8)	Not reported	32 (IQR 23-44)	Mean 48.8 (SD 14.4)	2.6 (IQR 1.1–6.8)	36 (IQR not reported)	
Total partici- pants	(	343	63	123	300	165	131	
Population	I from previous page	Adults aged ≥18 years. Some patients with presumptive tuberculosis, come with come with confirmed positive Gene Xpert sputum testing	Adults diagnosed with pulmonary tuberculosis	Adults aged ≥18 years with culture confirmed tuberculosis	Adults with a clinical diagnosis of pulmonary tuberculosis	Children aged <15 years with presumptive pulmonary tuberculosis	Adults aged ≥16 years with HIV and a known pulmonary tuberculosis diagnosis pre- treatment or symptoms of tuberculosis	
	(Continued	Luabeya 2019, South Africa <sup>19</sup>	Mesman 2019, Peru <sup>is</sup>	Mesman 2020, Peru <sup>17</sup>	Molina- Moya 2020, Moldova <sup>31</sup>	Nicol 2019, South Africa <sup>24</sup>	Shapiro 2022, Africa <sup>7</sup>	

e of NAAT		LAMP assay	nouse PCR plied systems oOnePlus Real- e PCR system)	iouse PCR plied systems oOnePlus Real- e PCR system)	
TYF		d TB- °C	it (Ap Bio: Stel	it (Ap Bio: Stel	
Laboratory processing		Sample centrifuged then resuspended. Pellet heated at 90° in extraction solution	DNA extracted usin QIAmp DNA miniki	DNA extracted usin QIAmp DNA miniki	
Sample storage		Not specified	Frozen and stored at -80°C	Frozen and stored at -80°C within 8 h of collection	
Sample transport		Sterile normal saline	In house lysis buffer	Sterile lysis buffer	temperature.
Number, location of swabs		3 tongue swabs (morning, spot, and evening)	3 buccal swabs	2 tongue swabs	ication. Temp≓
Swab type		Not specified	Whatman OmniSwab	Puritan Purflock Ultra and Copan FLOQSwabs	othermal amplii
Positive by oralswab and reference standard		38	18	44	op-mediated is
Reference standard (primary)		Sputum Xpert MTB/ RIF or culture	Sputum Xpert MTB/ RIF	Sputum Xpert MTB/ RIF Ultra or culture	=tuberculosis lc
Participants living with HIV (%)		0	0	36 (25%)	ogroup. TB-LAMP
Children aged <5 years		AN	AN	NA	epending on sub
Median age (years)		48-5 (range 17–88)	Median 38-45*	Median 30-34*	cation test. *D
Total partici- pants	(ab	101	40	194	ic acid amplifi
Population	d from previous pa	Adults aged >16 years with presumptive pulmonary tuberculosis	Adults aged ≥21 years with positive sputum Gene Xpert	Adults aged >18 years with presumptive pulmonary tuberculosis	licable. NAAT=nucle
	(Continue	Song 2021, China³²	Wood 2015, South Africa <sup>16</sup>	Wood 2021, Uganda <sup>18</sup>	NA=not app

Table: Characteristics of oral swabs for the diagnosis of pulmonary tuberculosis by study

Of the 20 studies included, 10 (50%) had low risk of bias in all four QUADAS-2 domains. Seven studies (35%) were considered high risk of bias for patient selection due to using a case-control design.<sup>14-19</sup> One study was considered high risk of bias in the index test domain due to incomplete description of oral swab collection and processing methods, and not stating a prespecified threshold for a positive index test.20 One study was considered high concern for applicability for patient selection, as it was carried out in a tertiary care referral hospital; however the study was carried out in a population of children with a relatively high prevalence of HIV (17.1%), both of which are target population characteristics for non-sputum-based diagnostic samples.<sup>21</sup> All studies were assessed as unclear concern for applicability of the index test as there is no standard methodology defined for using oral swabs for tuberculosis detection. All studies were considered low risk of bias in the reference test domain, and 18 were assessed as low risk of bias in the flow and timing domain.7.14-24,27,30-32 Two studies were assessed as unclear risk of bias as the interval between index and reference testing was not specified (see appendix pp 7–10 for further details).<sup>27</sup>

The 20 included studies from 16 reports were conducted in eight countries, including South Africa (seven studies), Peru (five), Uganda (three), South Korea (one), Kenya (one), Brazil (one), China (one), and Moldova (one). In adults, two (10%) studies were conducted in populations with a HIV prevalence of more than 50%.<sup>723</sup> One study used oral swabs as a screening specimen in a high-risk, asymptomatic population and all other studies evaluated oral swabs as a diagnostic specimen for participants with signs or symptoms of pulmonary tuberculosis.<sup>14</sup> Key elements of the included studies such as type of swab used, number of swabs collected per participant, site of the oral cavity swabbed, sample processing details, and type of NAAT used, are summarised in the table.

15 studies (12 reports, 2223 participants) evaluated oral swabs for tuberculosis detection in adults. Sensitivity of oral swabs and molecular test for the diagnosis of pulmonary tuberculosis in adults against a microbiological reference standard ranged from 36% to 91%, with specificity of oral swabs ranging from 66% to 100% (figure 3A, appendix p 11).

Sensitivity of tongue swabs in adults (eight studies) ranged from 52% to 91%; sensitivity of buccal swabs in adults (six studies) ranged from 36% to 90% (figure 3). The study with the lowest sensitivity for tongue swabs (52%) was performed in asymptomatic individuals as a screening test.<sup>14</sup> When used as a diagnostic test for individuals with symptoms of pulmonary tuberculosis, sensitivity of tongue swabs in adults ranged from 66% to 92% with specificity ranging from 66% to 100%.

Swab types evaluated included Whatman OmniSwab (six studies [30%]), Copan FLOQSwabs (five studies [25%]), Whatman EasiCollect FTA (two studies [10%]), Puritan PurFlock Ultra (two studies [10%]), multiple



Figure 3: Forest plot (A) and SROC plot (B) of oral swabs in adults using a microbiological reference standard by location of oral cavity swabbed OSA=oral swab analysis. SROC=summary receiver operating characteristic.

swab types (two studies [10%]), PrimeSwab Flocked Swab (one study [5%]), OmniGene Oral OMR-110 Kit (one study [5%]), and unspecified type (one study [5%]). Studies using Copan FLOQSwabs had sensitivities ranging from 66% to 91% and specificity from 66% to 100%; sensitivity using Whatman OmniSwabs ranged from 39% to 90% and specificity from 83% to 100% (appendix p 11).

NAATs used included manual PCR (14 studies [70%]), Xpert MTB/RIF Ultra (five [25%]), and tuberculosis loopmediated isothermal amplification (TB-LAMP; one [5%]). In-house PCR had sensitivities ranging from 36% to 91% and specificities from 66% to 100%, whereas studies using Xpert MTB/RIF Ultra had sensitivities ranging from 45% to 78% with all specificities reported as 100% (appendix p 11). In the most recent study using Xpert MTB/RIF Ultra, Andama and colleagues achieved a sensitivity of 78% (95% CI 64–88) and a specificity of 100% (95% CI 97–100) using tongue swabbing and new sample processing methods using Cepheid sample reagent.<sup>22</sup>

Four studies in adults reported invalid swab results (13 [1 $\cdot$ 8%] of 727 total swabs), including six in LaCourse and colleagues, four in Mesman and colleagues, two in Luabeya and colleagues, and one in Andama and colleagues.<sup>15,19,22,23</sup> Other studies either reported no invalid swab results or did not report the data on invalid results.

Two studies were conducted in populations with an HIV prevalence of greater than 50%.<sup>7,23</sup> Shapiro and



Figure 4: Forest plot of oral swabs in children using a composite (A) and microbiological (B) reference standard

colleagues (tongue swab, COPAN FLOQSwabs, in-house PCR) had an HIV prevalence of 92% and showed a sensitivity of 66% (95% CI 53–77) and a specificity of 78% (95% CI 66–87) using a quantification cycle (Cq) threshold of 38.<sup>7</sup> LaCourse and colleagues (buccal swab, Whatman OmniSwab, in-house PCR) had an HIV prevalence of 54%, a sensitivity of 68% (43–87), and a specificity of 83% (72–90; figure 3A).<sup>23</sup>

Five studies (four reports, 815 participants) evaluated sensitivity of oral swabs as a specimen in children. Sensitivity in children ranged from 8% to 42% when using a microbiological reference standard, and 5% to 42% when using a composite reference standard (figure 4). Specificities ranged from 93% to 100%. Only one study in children reported invalid results, with four invalid swab results.<sup>24</sup> Studies varied widely in terms of location of oral cavity swabbed, swab type, number of swabs collected, swab processing method, and molecular test used. All these factors impact sensitivity and specificity of oral swabs and molecular testing for the diagnosis of pulmonary tuberculosis. Given this, no meta-analysis of results was performed as studies differed too greatly.

# Discussion

In this systematic review of the diagnostic accuracy of oral swabs and molecular testing for pulmonary tuberculosis, we found that the sensitivity of oral swabs ranged from 36% to 91% for adults and between 5% and 42% for children. As expected, sensitivity was reduced in children who are more likely to have paucibacillary disease and an imperfect reference standard. Sensitivity of oral swabs as a clinical specimen varies in both adults and children when diverse methods are used. Specificity is less variable than sensitivity, with most studies reporting specificity greater than 90%.

Optimising sample processing methods, including swab type, anatomical location sampled, and NAAT platforms could lead to increased diagnostic accuracy with oral swabs. Luabeya and colleagues demonstrated that Cq values for MTB DNA were significantly lower (corresponding to stronger PCR signals) in samples collected from the tongue compared with the cheek or gingiva.<sup>22</sup> Copan FLOQSwabs showed increased biomass collection relative to Whatman OmniSwabs and Puritan Purflock Ultra Swabs, which has been linked to higher sensitivity.<sup>19</sup> Andama and colleagues in Uganda evaluated a protocol for tongue swabbing coupled with Xpert MTB/ RIF Ultra using new sample processing methods. This was the first study in adults with presumptive tuberculosis to demonstrate high sensitivity of oral swabs using Xpert MTB/RIF Ultra. The authors found overall sensitivity of 72% (95% CI 59-83) and specificity of 100% (95% CI 97-100) using a microbiological reference standard,<sup>21</sup> achieved using a double swabbing strategy and processing with Cepheid sample reagent and Xpert MTB/RIF Ultra. Using an alternative processing method, where a single swab was boiled, incubated, and mixed without using sample reagent, sensitivity was very similar at 77% (60–90) and specificity was 100% (16-100).21 With the increased biomass collected by double swabs, the processing method using Cepheid sample reagent showed similar sensitivity to the boil method (73% using the Cepheid sample reagent vs 77% using the boil method). The use of sample reagent and similar processing methods to sputum potentially allows laboratories already performing sputum testing with Xpert MTB/RIF Ultra to use oral swabs as a sample type with minimal change in processing methods.

Several other differences between studies included uses of different media types for swab transport, freezing of swab samples before testing versus immediate testing, and laboratory processing before PCR testing. These methods are continuously being evaluated but are not yet optimised. A consensus platform protocol for future research studies of oral swabs, which standardises procedures for swab collection, processing, storage, and other procedures, would improve consistency within studies, facilitate comparisons between studies, and allow more accurate accuracy estimates.

We were unable to perform meta-analyses as originally intended in the protocol because of considerable heterogeneity related to location of sampling, swab type, transport and storage methods, and type of NAAT used, as well as a limited number of studies for each combination of these factors (figure 1). A recent study meta-analysed several of the studies included in this review, and showed a summary sensitivity of 58.5% and specificity of 85.6% for oral swabs in the diagnosis of pulmonary tuberculosis.<sup>25</sup> However, the sources of heterogeneity we identified were not investigated as part of that analysis.

Most studies included in this review had low risk of bias in all four QUADAS-2 domains. However, methods for evaluating accuracy of oral swabs were inconsistently described, which means some studies lack details of sample collection, transport, storage, and processing that are important to reproduce the results without substantial heterogeneity. This could lead to unrecognised sources of variability in sensitivity and specificity. In particular, there was frequent use of in-house PCR assays in many of the studies, which can vary considerably in specificity when compared with commercial PCR assays.<sup>26</sup> In addition, many of these studies were conducted in populations with a high pre-test probability of disease and would be expected to have higher bacillary burdens, which is associated with increased sensitivity. Two studies showed that smear status and cavitary disease are associated with higher sensitivity.<sup>17,27</sup> Regarding applicability of our findings to the review question, in the index test domain, we considered applicability to be unclear as, currently, there is not a standardised protocol for oral swabs. The lack of a standardised protocol, with each study using a different approach, limited our ability to assess diagnostic accuracy by the methods evaluated in the studies.

The findings in this review were based on comprehensive searching, strict selection criteria, and standardised data extraction. To identify studies, we searched multiple databases up to July 5, 2022, without language restriction. However, we could have missed studies despite the comprehensive search. We corresponded with primary study authors to obtain additional data and information that was missing from the papers. In addition, this is a rapidly developing field, and a modified search update identified additional evidence published or submitted since the full search was completed on Aug 23, 2023. Although no additional studies using swabs have been published, additional studies<sup>28</sup> have continued to explore saliva with findings of good sensitivity comparable with the highest sensitivities seen with oral swabs.

For detection of pulmonary tuberculosis, most studies evaluated oral swab specimens collected from participants with presumptive tuberculosis attending primary care facilities and local hospitals. Hence, for most studies, the participant characteristics and settings matched our review question. One study evaluated oral swab specimens as a screening test for pulmonary tuberculosis in an asymptomatic population at high risk for infection.<sup>14</sup> Although this is a potential use of oral swabs for tuberculosis detection, it did not reflect the primary purpose of this review; further study in community-based settings and asymptomatic populations are needed to clarify these roles. Of note, since the majority of studies in adults required sputum to be collected for a microbiological reference standard, individuals unable to provide a sputum sample might have been underrepresented in this review. As sputum-scarce populations could particularly benefit from oral swab sampling, this is a limitation in the current literature on oral swabs which might misestimate the diagnostic yield of oral swabs in important target populations. As higher bacillary load does seem to predict increased sensitivity, it is crucial to optimise methods for sputum-scarce populations such as children and people living with HIV as they frequently present with lower pulmonary bacillary loads.

We found that sensitivity of oral swabs varies in both adults and children when diverse methods are used. Specificity is less variable than sensitivity with most studies reporting specificity greater than 90%. As oral swabs are not currently approved or recommended for use in the diagnosis of pulmonary tuberculosis, this review does not alter current clinical practice. We recommend that research going forward compares specific aspects of oral swabs, either through prospective cohort or clinical trials, so that optimal methods can be determined, such as direct comparisons of two swab types or of two processing methods such as seen in Andama and colleagues.<sup>22</sup>

Oral swabs with molecular testing can provide accurate results for the diagnosis of pulmonary tuberculosis. The highest sensitivities in this review approach the acceptable sensitivity range (as defined by WHO consolidated guidelines on tuberculosis<sup>2</sup>) for an initial diagnostic test relative to performance of sputum-based diagnostics, suggesting that with continued optimisation, swabs could be a good sample type for tuberculosis detection.29 Future studies should use a standardised protocol for comparing different aspects of oral swab collection and processing, keeping all aspects of the protocol the same except one. To advance the potential of oral swabs in a point-of-care diagnostic test, further study using existing point-of-care molecular WHO-recommended rapid diagnostics should be considered, as well as development of novel, lowcomplexity point-of-care platforms for use with oral swabs.

## Contributors

ECC, AES, and KRS contributed to study design, literature search and review, data extraction, data analysis, figure creation, and manuscript writing. GAC and MK contributed to study design, data analysis, figure

creation, and manuscript writing. MR contributed to study design and manuscript writing. All authors had access to the data for this study. ECC and AES independently accessed and verified all data. All authors have read and approved the final version of this manuscript and had final responsibility for submission of the manuscript.

### **Declaration of interests**

This study received funding from FIND. ECC reports institution payments from FIND. MK and MR are employed by FIND. AES is supported in part by a National Institutes of Health (NIH) K23 AI40918 award. GAC is funded by the NIH and Bill and Melinda Gates foundation and reports receiving donations of research supplies (FLOQswabs) from Copan Italia. KS has received financial support from Cochrane Infectious Diseases, McGill University, Baylor College of Medicine, Maastricht University, TB Proof, and WHO Global Tuberculosis Programme; consultancy fees from FIND, the global alliance for diagnostics; consulting fees from Stellenbosch University, and travel support to attend WHO guideline development group meetings.

#### Data sharing

The data included in this review are all published data or available by pre-print through authorised subscriptions. The protocol for this study is not available on any protocol storage site but is available on request from the corresponding author.

#### Acknowledgments

We are grateful to Ethan Valinetz of Indiana University School of Medicine, Indianapolis, IN, USA for content input on oral swabs for infectious diseases testing and Vittoria Lutje of Cochrane Infectious Diseases Group, Liverpool, UK for informatics support for this study.

#### References

- 1 WHO. Global tuberculosis report 2023. Geneva: World Health Organization, 2023.
- 2 WHO. WHO consolidated guidelines on tuberculosis. Module 3: diagnosis. Rapid diagnostics for tuberculosis detection, 2021 update. https://www.who.int/publications/i/item/9789240030589 (accessed July 5, 2022).
- 3 Zifodya JS, Kreniske JS, Schiller I, et al. Xpert Ultra versus Xpert MTB/RIF for pulmonary tuberculosis and rifampicin resistance in adults with presumptive pulmonary tuberculosis. *Cochrane Database Syst Rev* 2021; 2: CD009593.
- 4 Kay AW, González Fernández L, Takwoingi Y, et al. Xpert MTB/RIF and Xpert MTB/RIF Ultra assays for active tuberculosis and rifampicin resistance in children. *Cochrane Database Syst Rev* 2020; 8: CD013359.
- 5 Valinetz ED, Cangelosi GA. A look inside: oral sampling for detection of non-oral infectious diseases. *J Clin Microbiol* 2021; 59: e02360-20.
- 6 Codsi R, Errett NA, Luabeya AK, et al. Preferences of healthcare workers using tongue swabs for tuberculosis diagnosis during COVID-19. PLoS Glob Public Health 2023; 3: e0001430–0.
- 7 Shapiro AE, Olson AM, Kidoguchi L, et al. Complementary nonsputum diagnostic testing for tuberculosis in people with HIV using oral swab PCR and urine lipoarabinomannan detection. *J Clin Microbiol* 2022; 60: e0043122.
- 8 WHO. Definitions and reporting framework for tuberculosis 2013 revision: updated December 2014 and January 2020. April 9, 2013. https://www.who.int/publications/i/item/9789241505345 (accessed June 15, 2022).
- 9 WHO. WHO operational handbook on tuberculosis. Module 5: management of tuberculosis in children and adolescents. March 16, 2022. https://www.who.int/publications/i/item/ 9789240046832 (accessed June 15, 2022).
- 10 Covidence systematic review software [Computer program]. Covidence. Melbourne, VIC: Veritas Health Innovation. 2015. https://www.covidence.org/ (accessed Jan 27, 2020).
- 11 Deeks JJ, Bossuyt PM, Leeflang MM, Takwoingi Y, eds. Cochrane Handbook for systematic reviews of diagnostic test accuracy. Version 2.0 (updated July 2023). Cochrane, 2023. https://training.cochrane.org/handbook-diagnostic-test-accuracy/ current (accessed Aug 23, 2023).

- 12 Whiting PF, Rutjes AW, Westwood ME, et al. QUADAS-2: a revised tool for the quality assessment of diagnostic accuracy studies. *Ann Intern Med* 2011; 155: 529–36.
- 13 Homorodean D. Compararea rezultatelor culturilor obținute în mediu lichid în sistemul MB/BacT şi pe mediul Löwenstein Jensen de la pacienți selectați. *Pneumologia* 2005; 54: 62–65 (in Romanian).
- 14 Lima F, Santos AS, Oliveira RD, et al. Oral swab testing by Xpert® MTB/RIF Ultra for mass tuberculosis screening in prisons. *J Clin Tuberc Other Mycobact Dis* 2020; **19**: 100148.
- 5 Mesman AW, Calderon R, Soto M, et al. Mycobacterium tuberculosis detection from oral swabs with Xpert MTB/RIF ULTRA: a pilot study. BMC Res Notes 2019; 12: 349.
- 16 Wood RC, Luabeya AK, Weigel KM, et al. Detection of *Mycobacterium tuberculosis* DNA on the oral mucosa of tuberculosis patients. *Sci Rep* 2015; 5: 8668.
- 17 Mesman AW, Calderon RI, Pollock NR, et al. Molecular detection of Mycobacterium tuberculosis from buccal swabs among adult in Peru. Sci Rep 2020; 10: 22231.
- 18 Wood RC, Andama A, Hermansky G, et al. Characterization of oral swab samples for diagnosis of pulmonary tuberculosis. *PLoS One* 2021; 16: e0251422.
- Luabeya AK, Wood RC, Shenje J, et al. Noninvasive detection of tuberculosis by oral swab analysis. *J Clin Microbiol* 2019; 57: e01847–18.
- 20 Kang YA, Koo B, Kim OH, et al. Gene-based diagnosis of tuberculosis from oral swabs with a new generation pathogen enrichment technique. *Microbiol Spectr* 2022; 10: e0020722.
- 21 Ealand C, Peters J, Jacobs O, et al. Detection of *Mycobacterium tuberculosis* complex bacilli and nucleic acids from tongue swabs in young, hospitalized children. *Front Cell Infect Microbiol* 2021; 11: 696379.
- 22 Andama A, Whitman GR, Crowder R, et al. Accuracy of tongue swab testing using Xpert MTB-RIF Ultra for tuberculosis diagnosis. *J Clin Microbiol* 2022; 60: e0042122.
- 23 LaCourse SM, Seko E, Wood R, et al. Diagnostic performance of oral swabs for non-sputum based TB diagnosis in a TB/HIV endemic setting. *PLoS One* 2022; 17: e0262123–3.
- 24 Nicol MP, Wood RC, Workman L, et al. Microbiological diagnosis of pulmonary tuberculosis in children by oral swab polymerase chain reaction. *Sci Rep* 2019; **9**: 10789.
- 25 Savage HR, Rickman HM, Burke RM, et al. Accuracy of upper respiratory tract samples to diagnose *Mycobacterium tuberculosis*: a systematic review and meta-analysis. *Lancet Microbe* 2023; 4: e811–21.
- 26 Wei Z, Zhang X, Wei C, et al. Diagnostic accuracy of in-house realtime PCR assay for *Mycobacterium tuberculosis*: a systematic review and meta-analysis. *BMC Infect Dis* 2019; 19: 701.
- 27 Flores JA, Calderón R, Mesman AW, et al. Detection of Mycobacterium tuberculosis DNA in buccal swab samples from children in Lima, Peru. Pediatr Infect Dis J 2020; 39: e376–80.
- 28 Byanyima P, Kaswabuli S, Musisi E, et al. Feasibility and sensitivity of saliva GeneXpert MTB/RIF Ultra for tuberculosis diagnosis in adults in Uganda. *Microbiol Spectr* 2022; 10: e0086022.
- 29 Drain PK, Gardiner J, Hannah H, et al. Guidance for studies evaluating the accuracy of biomarker-based nonsputum tests to diagnose tuberculosis. J Infect Dis 2019; 220 (suppl 3): S108–15.
- 30 Cox H, Workman L, Bateman L, et al. Oral swab specimens tested with Xpert MTB/RIF ultra assay for diagnosis of pulmonary tuberculosis in children: a diagnostic accuracy study. *Clin Infect Dis* 2022; 75: 2145–52.
- 31 Molina-Moya B, Ciobanu N, Hernández M, et al. Molecular detection of *Mycobacterium tuberculosis* in oral mucosa from patients with presumptive tuberculosis. J Clin Med 2020; 9: 4124.
- 32 Song Y, Ma Y, Liu R, et al. Diagnostic yield of oral swab testing by TB-LAMP for diagnosis of pulmonary tuberculosis. *Infect Drug Resist* 2021; 14: 89–95.